

USER'S PERSONAL DATA****

Name

Surname

Phone

E-mail

Age

Pathology

Notes

WHEELCHAIR FEATURES

Brand

Model

Front tubes diameter

Seat width

Seat depth (b)

Toe to front tube distance (a)*

Ground/knee hight (d)

Lenght of stretched arms from the shoulder back (c)

Fron tubes inclination

User's weight


Footrest

Fixed

Internal

Extractable

External



For wheelchairs with removable footrests, measure as per "a" from toe to front tube.

